

tsRTarget

Requirements:

- Linux system, enough disk space and Ram depending on the size of RNA deep sequencing data. (Tested system: ubuntu 12.04 LTS, ubuntu 16.04 LTS)

Installation

- Download tsRTarget pipeline package from <https://github.com/zhlingl/tsRFun>

```
wget https://github.com/zhlingl/tsRFun/archive/refs/heads/main.zip
```

- Download necessary software, packages and reference databases as listed below:
 - The CLIP Tool Kit (CTK) (https://zhanglab.c2b2.columbia.edu/index.php/CTK_Documentation) Tested version:1.1.3)
 - RNAhybrid (<https://bibiserv.cebitec.uni-bielefeld.de/rnahybrid/>)
 - BLAST(<https://blast.ncbi.nlm.nih.gov/Blast.cgi>) (Tested version:2.10.1)
 - Reference database (See lists and download link of all pre-compiled species' databases in Pre-compiled Databases Instruction)
- Download tsRTarget pipeline package

- Install tsRTarget from <https://github.com/zhlingl/tsRFun>

- Unpack tsRTarget package

```
unzip main.zip
```

- Attach the tsRTarget directory to your PATH:

```
echo 'export PATH=$PATH:your_path_to_tsRFun-master/tsRTarget' >> ~/.bashrc.  
chmod 755 your_path_to_tsRTarget-master/tsRTarget/tsRTarget.sh
```

- Install CTK

- Unpack czplib-1.0.x.zip

```
unzip czplib-1.0.x.zip  
mv czplib-1.0.x /usr/local/lib/czplib
```

- Add the library path to the environment variable, so perl can find it.

```
export PERL5LIB=your_path_to_czplib
```

- Download CTK code and likewise decompress and move to whatever directory you like (as an example, we use /usr/local/)

```
unzip ctk-1.0.x.zip  
mv ctk-1.0.x /usr/local/CTK
```

- Attach the CTK directory to your PATH:

```
echo 'export PATH=$PATH:your_path_to_CTK'
```

- Install RNAhybrid

- Unpack RNAhybrid-2.1.2.tar.gz.

```
wget https://bibiserv.cebitec.uni-bielefeld.de/applications/rnahybrid/resources/downloads/RNAhybrid-2.1.2.tar.gz  
tar -xzf RNAhybrid-2.1.2.tar.gz  
cd /path/to/ RNAhybrid-2.1.2  
./configure  
make
```

```
make install
cd bin/
```

- Attach the bedtools directory to your PATH:

```
echo 'export PATH=$PATH:your_path_to_RNAhybrid' >> ~/.bashrc
```

- Start a new shell session to apply changes to environment variables:

```
source ~/.bashrc
```

- Install BLAST

- Unpack ncbi-blast-2.10.1+-x64-linux.tar.gz from <https://ftp.ncbi.nlm.nih.gov/blast/executables/blast+/LATEST/>.

```
tar -zxvf ncbi-blast-2.10.1+-x64-linux.tar.gz
```

- Attach the samtools directory to your PATH:

```
echo 'export PATH=$PATH:your_path_to_BLAST' >> ~/.bashrc
```

- Start a new shell session to apply changes to environment variables:

```
source ~/.bashrc
```

- Download tsRFun pipeline package

- Test if everything is installed properly:

```
perl -v
bash tsRTarget.sh -h
ctk -h
RNAhybrid -h
blastn -h
```

- If you get any error messages you should install the software or perl modules once again.

Script description

- The input files of tsRTarget.sh are:

```
Options:
-i <inputfile>      Input could be:
                    a .fastq/.fq or .fasta/.fa file.
-t data type clear/clash or CLIP (iclip/eclip/par-clip/hits-clip)
-o outputdir        address of annotation results
-a index_address
-m mismatch         default=1
-w min tsRNA_length word_size default=14
-n max tsRNA_length word_size default=40
-g min target_length default=10
-h max -g min target_length default(70 for clear clash data; 140 for clip data)
-s matched-length   default=6
-e gapopen          default=1
-P Collapse PCR duplicates default=T
-S seed T/N         default=T
-R seed start       default=2
-E seed end         default=8
-M Min Free Energy  default=-10
-c match_clip       the word_size of RNA-RNA target; only used for CLIP seq data,

Others:
-h                  print this usage message
```

Example of use:

```
bash tsRTarget.sh -i example.fa -t clear -o example_result -a hg38_index
```